

4.2 SmartMan Manual – 4.Trainer Menu 2. Intubation Trainer

v5.1 and later Current version v5.3.3.3

Click on any row in **Table of Contents** to go to that section

Table of Contents

4.2 SmartMan Manual – 4.Trainer Menu 2. Intubation Trainer	1
4.0 The Trainer Menu	1
4. 2. Intubation Trainer	1
4.2.1 What it does	1
4.2.2 What you are to do	1
4.2.3 How to Start The Activity	2

4.0 The Trainer Menu

The Trainer Menu provides special activities for skills which need to be mastered in order to be able to perform CPR at a high level.

This area is undergoing development and there are only two skills at this time. As new activities are developed to improve specific components skills they will appear in the Trainer Menu

4.1 Ventilations Trainer

Here you practice improving how you perform a ventilation as required for rescue breaths and in CPR during 30:2. It displays volume as you enter air and a timer is running so that you learn how to squeeze about 600ml of air in 1 second.

It also shows the tidal nature of a breath so that you can make sure you fully release the BVM

4.2 Intubation Trainer

This is a training activity to help you improve the speed and accuracy with which you place an intubation device.

[Top](#)

4. 2. Intubation Trainer

4.2.1 What it does

This activity is to help improve how quickly you can establish an effective advanced airway. The airway manikin will accept most devices that are placed through the mouth into the airway.

This activity will time how long it takes for you to place the intubation device and verify that the cuff is properly seated and properly inflated.

[Top](#)

4.2.2 What you are to do



The protocol for this activity is to prepare the devices, start the activity, place the device, provide about 600ml of air and hold the BVM for 5 seconds. During this time the program will verify the seal.

At completion of the activity, the program will provide you a score on the effectiveness of placing and of the time it took.

[Top](#)

4.2.3 How to Start The Activity

Click through the selections with the mouse or by pressing the Space bar.



and



When you press “start” the timer is running. You will see skinny gray bars on the screen indicating that it is recording.

Place the intubation device as quickly as possible.


In the picture below

Gray Bars: The air started going into the lungs at 10 seconds.

Green: Air was put in to about 570ml

Top of Green bar: This is perfectly horizontal. The cuff held perfectly.

The summary of the results is in the Session Results area on top.



Press Re-try to practice another intubation.

Press Menu to return to the Main Menu.

[Top](#)